

# META-ANALYSIS: COMBINING NPWT + TOT

11 RCTS  
844 PATIENTS

Su J, Zhang D, Du J, et al. The Efficacy of Negative Pressure Wound Therapy Combined with Topical Oxygen Therapy in Treating Chronic Refractory Wounds: A Systematic Review and Meta-Analysis. *Adv Wound Care*. Published online August 13, 2025.

## Study Objective

To evaluate the clinical efficacy of combining **Negative Pressure Wound Therapy (NPWT)** with **Topical Oxygen Therapy (TOT)** for treating chronic refractory wounds (CRWs).

**Rationale:** NPWT often induces local hypoxia. TOT supplements oxygen to restore metabolic function and optimize the wound microenvironment.

1.51×

Healing Rate  
(RR: 1.51)

73%

Lower Infection  
(RR: 0.27)

-9.1

Days to Heal  
(Mean Diff)

## Methodology

- Design:** 11 Randomized Controlled Trials (RCTs)
- Participants:** 844 patients with chronic refractory wounds
- Protocol:** PRISMA 2015 guidelines (PROSPERO: CRD42024618561)
- Data Sources:** PubMed, Cochrane, Embase, Web of Science, CNKI (through Oct 2024)

## Meta-Analysis Results Summary

Comparison: NPWT+TOT vs. NPWT Alone

Complete Wound Healing Rate	RR	1.51	1.36 – 1.69
Bacterial Positivity Rate	RR	0.27	0.18 – 0.41
Overall Healing Time	MD (days)	-9.09	-11.98 to -6.20
Time to Skin Grafting	MD (days)	-2.82	-3.15 to -2.50
Granulation Coverage	MD (%)	+7.56	6.09 to 9.03

RR = Risk Ratio; MD = Mean Difference; CI = Confidence Interval

## Included Wound Etiologies:

Diabetic Foot Ulcers

Pressure Injuries

Venous Ulcers

Mixed Etiologies

# Key Outcomes & Conclusion

## Subgroup Analysis & Clinical Implications

### Subgroup Analysis: Healing Rates

Pressure Injuries



RR 1.89 ✓

Highest Efficacy

Mixed CRWs



RR 1.38 ✓

Significant Benefit

Diabetic Foot



RR 1.75 (Trend Only)

Small Sample Size

### Why NPWT + TOT Works



#### Reverses Hypoxia

Counteracts oxygen deficiency from negative pressure suction.



#### Promotes Angiogenesis

Stimulates VEGF expression to build new capillaries.



#### Bactericidal Effect

Enhances neutrophil respiratory burst to clear infection.



#### Accerelates Repair

Supports mitochondrial energy for tissue regeneration.

### Clinical Benefits

- 🕒 Significant reduction in time to skin grafting (-2.82 days)
- 🧼 Cleaner wound bed with significantly lower bacterial load
- 💰 Potential for reduced hospital stay and resource utilization
- 🔄 Helps transition stuck wounds from inflammatory to proliferative phase



### Conclusion

Adding Topical Oxygen Therapy to NPWT is a synergistic strategy that significantly improves healing rates and reduces infection risks in chronic refractory wounds, particularly pressure injuries.

### Recommendation

Consider adopting combined protocol for wounds that fail to progress with standard NPWT alone.